


ICWUSA.com Product Bulletin




The purpose of the ICW Product Bulletin is to share interesting solutions using ICW products. Watch for new bulletins by issue number and date. Past issues are available on our website: www.icwusa.com/product-bulletins/

Service Trays

SRVCTRY



KHT



ICW service trays can mount on all of ICW's mounting configurations. We have five different sizes of trays that can hold anything from pens to printers. Additionally, ICW can custom build trays for those looking for specific dimensions.


Our trays can be mounted to ICW mounts, wall tracks or directly to the wall. Most trays can swivel horizontally and flip up to stow out of the way. We also have keyboard work surfaces that double as tray mounts.

Available Service Tray Sizes:

- TRY1611 (16.75" x 11.75")
- TRY1718 (17.6" x 18.3")
- SCANTRAY (6" x 8")
- SRVCTRY (16.2" x 5.7")
- KHT (10.2" x 12.2")

For more information on ICW solutions, please email sales@icwusa.com or call 1-800-558-4435 to speak to one of our product specialists.

SCANTRAY



TRY1611

