

Elite Laptop Arm



Model Number as shown:
EL5220LT-WT-QG
(laptop not included)

ICW

Elite Laptop Arm

DESCRIPTION:

The double arm Elite with a laptop tray incorporates design features allowing it to support up to 28 lbs (12.7 kg) while providing effortless movement. It allows a laptop to be pulled out when required. The Elite arm is fully cable managed, gas spring assisted, and provides up to 42" of vertical and horizontal adjustment from mount to tray, along with 180° of rotation, making it ideal in medical operating rooms or dental operatories.


Machined from billet aluminum, the unit is powder coated for a durable, lasting finish. Hand assembly ensures that a quality product is shipped to you.

FEATURES:


- Fully enclosed gas cylinders
- Hardened alloy shoulder bolts and Teflon bushings, engineered for years of reliable service and quiet operation
- Powder coated and anodized for long wear and easy cleaning
- Cable channeling through the arm assembly
- Flip-up tray

SPECIFICATIONS:

Finish:	Gray, black or medical white
Weight capacity:	28 lbs (12.7 kg)
Vert. adjustment range:	42"
Maximum reach:	66"
Arm swivel at mount:	180°
Tray swivel:	360°
Tray size:	12" x 15"




WALL PLATE




SIDE

FOLDED TO
MINIMUM PROFILE



FRONT



SIDE, EXTENDED

